





OverdenSURE®

Attachment System

Technique Manual



 **ZimVie**



OverdenSURE

Attachment System

Get connected with ZimVie's removable attachment system, OverdenSURE!

ZimVie offers an extensive OverdenSURE Removable Attachment System compatible with our ZimVie Implant solutions plus 12 of the industry's leading systems! Whatever you need for removable full-arch restoration cases, we've got your connection!

ZimVie Implant Systems

- TSX® Implant System
- TSV® and Eztetic® Implant Systems
- Certain® Implant Systems
- External Hex® Implant System
- Tapered Swiss Plus® Implant System

Compatible with the following Systems

- Ankylos® C/X
- Astra Tech® EV/Osseospeed® TX
- BioHorizons® Internal
- Camlog®
- NobelActive® /NobelReplace®
- Nobel/Replace Select™ Tapered
- Nobel Brånemark
- Straumann® Bone/Tissue level
- Neodent® Grand Morse™



Table of Contents

System Overview – OverdenSURE Abutments	4
General Information	4
Selecting the OverdenSURE Abutments	4
Selecting the OverdenSURE Retention Inserts	4
Specifications	4
Instrumentation and Components Needed	5
Seating the OverdenSURE Abutments	6
Tightening the OverdenSURE Abutments	6
Confirm the OverdenSURE Abutments are Seated	6
Direct Restorative Technique	7
Prepare the Denture for Processing	7
Pick-Up OverdenSURE Housing Caps	8
Remove OverdenSURE Black Inserts	8
Place Final OverdenSURE Retention Inserts	8
Deliver the Removable Prosthesis	9
Indirect Restorative Technique	10
Clinician – Take Impressions of the OverdenSURE Abutments	10
Transfer Technique	11
Laboratory – Fabricate the Working Cast	11
Laboratory – Record Base and Occlusion Rim	11
Clinician – Occlusion Record	11
Laboratory – Articulate	11
Clinician – Denture Try-In	12
Laboratory – Process the Removable Prosthesis	12
Clinician – Remove OverdenSURE Black Inserts	12
Clinician – Place Final OverdenSURE Retention Inserts	13
Deliver the Removable Prosthesis	13
Ordering Info	14


Note: Images shown throughout this manual are representational in nature and may not be to scale or display the exact geometry or color of the actual components.

OverdenSURE System Overview

OverdenSURE Abutments

Considerations

- Tissue-supported removable overdentures on two or more implants.
- Partially edentulous overdentures with one or more implants.
- Limited interarch distance.
- Angle correction with up to 20° between divergent implants when using the standard retention inserts (blue, pink, and clear) or the standard retention insert kits.
- Angle correction with up to 40° between divergent implants when using the extended retention inserts (red, yellow, and green) or the extended retention insert kits.



Selecting the OverdenSURE Abutments

- Choose a cuff length based on the tissue depth.
- Determine the tissue depth by measuring from the coronal rim of the implant body to the crest of the soft-tissue in three or four areas.
- Choose the corresponding abutment cuff height that equals the tissue measurement or is the next closest higher size.









Selecting the OverdenSURE Retention Inserts

The OverdenSURE Retention Inserts have been designed to provide the correct amount of strength for each level of retention while maintaining a color coding system that is intuitive to the clinician.

Specifications

- Abutment: Titanium alloy with light-gold Zirconium Nitride coating
- Housing: Titanium alloy
- Retention Inserts: Nylon

Retention Inserts			
Up to 20° Angle Correction		Up to 40° Angle Correction	
	1.5lbs / 680g		1.0lbs / 450g
	3.0lbs / 1360g		2.0lbs / 910g
	5.0lbs / 2270g		4.0lbs / 1810g

Instrumentation and Components Needed

- OverdenSURE Abutment Insertion and Removal Tool [ODS-IRTOOL]



- ZimVie Restorative Torque Wrench [TWR]



- OverdenSURE Abutment Driver [ODS-DRVR]



- OverdenSURE Impression Copings [ODS-IC01]



Seating the OverdenSURE Abutments

- Seat the abutment into the implant and thread it until finger-tight using the transporter included with each abutment.

Tightening the OverdenSURE Abutments

- Manually engage the OverdenSURE driver [ODS-DRVR] with the abutment, aligning the axis of the driver with the axis of the abutment and fitting the fingers of the driver over the abutment's circumference.


Tip: For easier driver-to-abutment engagement, slightly rotate the driver over the abutment until you feel the internal tri-lobe engage first, then apply normal pressure downward to ensure the driver's outer fingers have also engaged the abutment's circumference

- Insert the Torque Ratchet Wrench [TWR] onto the driver and torque the abutments to the indicated force:
 - Certain and ExHex OverdenSURE Abutments:
20 Ncm
 - TSX, Eztetic, and Tapered Screw-Vent® OverdenSURE Abutments:
30 Ncm

Caution: To remove the driver from the abutment, jiggle the driver slightly from side to side in order to disengage the driver from the abutment first. Taking care that the wrench is properly engaged with the driver, then remove the wrench [TWR] and driver [ODS-DRVR] from the mouth.

Confirm the OverdenSURE Abutments Are Seated

- Radiograph the interfaces to verify complete seating of the abutments on all of the implants.
- Place the film perpendicular to the interface of the abutment on the implant.



Direct Restorative Technique

Prepare the Denture for Processing

- Place the OverdenSURE Housing Caps [ODS-RH01] with the pre-loaded Black Inserts on the abutments and press down firmly to make sure they are fully seated.
- Place a Block-out Spacer [ODS-BS020] at the base of the abutments to prevent acrylic resin from flowing into undercuts of the Housing Caps.




Note: In all geographies with the exception of the U.S., an OverdenSURE Standard Abutment Retention Kit [ODS-RIKIT2001] is included with every abutment. One Housing Cap [ODS-RH01] pre-loaded with a Black Insert and one Block-out Spacer are included within this kit and packaged with every abutment.

Note: In the U.S., OverdenSURE Abutments are packaged individually and retention inserts (ODSZ-RIKIT2001 or ODSZ-RIKIT4001) must be ordered separately.

- Apply contact-detecting material to the intaglio surface of the denture. Seat the denture over the Housing Caps [ODS-RH01] and verify that it is completely seated on the ridge.



- With a recess bur create the relief pockets where the Housing Caps will fit in order to be processed into the denture.
- Create undercuts around the relieved pockets to increase mechanical retention and open lingual/palatal vents to allow excess material to flow out of the denture.




Pick-Up OverdenSURE Housing Caps

- Place autopolymerizing acrylic resin or light cure composite resin over the OverdenSURE Housing Caps [ODS-RH01] and into the relieved pockets in the denture.
- Place and position the denture in the mouth and have the patient close lightly into centric occlusion. Allow the acrylic resin to set per the manufacturer's instructions.
- Detach the denture from the abutments and verify that the Housing Caps have been successfully processed into the denture.
- Fill any voids with autopolymerizing acrylic resin or light cure composite resin. Once cured, remove any excess material with a grind bur.




Remove OverdenSURE Black Inserts

- Using the OverdenSURE Removal and Insertion Tool [ODS-IRTOOL], remove the Black Inserts from each Housing Cap.
- To do so, make sure the Removal and Insertion Tool is set to "removal mode."
- This is achieved by turning the golden head of the tool two full rotations counter-clockwise. You will notice the blue tip recess into the tool.
- Insert the tool into the Black Insert and tilt it so the edge of the tip engages the inner lip of the insert.
- Once the tool has full engagement, pull it out.
- To remove the Black Insert from the tool, turn the golden head of the tool two full rotations clockwise. The blue tip will emerge and push off the Black Insert.



Place Final OverdenSURE Retention Inserts

- Ensure the Removal and Insertion Tool [ODS-IRTOOL] is set to "insertion mode." To do so, completely unscrew the golden head from the blue portion of the tool.
- Select the appropriate Retention Insert and place it on the insertion tip of the tool. Firmly insert each final Retention Insert into the empty Housing Caps and verify they sit flush with the rim of the Housing Caps.



Note: In all geographies with the exception of the U.S., an OverdenSURE Standard Abutment Retention Kit [ODS-RIKIT2001] is included with every abutment.

One Blue Retention Insert (Extra Light), one Pink Retention Insert (Light), one Clear Retention Insert (Medium), and one Red Retention Insert (Extra Light; Extended Range) are included within this kit and packaged with every abutment.

Note: In the U.S., OverdenSURE Abutments are packaged individually and retention inserts (ODSZ-RIKIT2001 or ODSZ-RIKIT4001) must be ordered separately.

Deliver the Removable Prosthesis

- Seat the prosthesis on the OverdenSURE Abutments and press down so that the Retention Inserts snap onto the abutments.
- Verify the occlusion.
- Have the patient remove and insert the prosthesis to ensure the retention is adequate. If the retention level needs to be adjusted, select the appropriate Retention Insert and repeat the necessary steps.



Indirect Restorative Technique

Clinician – Take Impressions of the OverdenSURE Abutments

- A custom or stock closed impression tray may be used for this. Provide relief for the height of the OverdenSURE Abutments and Impression Copings [ODS-IC01] plus 2 mm. The Impression Copings are 4.5 mm in height.
- Place the OverdenSURE Impression Copings onto the abutments. Inject light-body impression material around the copings and record a full-arch impression with medium-body material. Allow the material to set per the manufacturer's instructions.
- Sometimes it is helpful to place impression adhesive on the impression copings to ensure their removal in the impression.

Note: *OverdenSURE Impression Copings [ODS-IC01] are sold separately.*

- Remove the impression from the mouth and verify that the impression material completely adapted around each OverdenSURE Impression Coping. The impression copings should remain inside the impression.
- An optional technique is to use the patient's existing denture in combination with the OverdenSURE Housing Caps [ODS-RH01] with a pre-loaded Black Insert for a reline impression.



Transfer Technique

- Insert the OverdenSURE Abutment Analogs [ODS-AA01] into the Impression Copings located within the impression material and send the impression to the laboratory.



Laboratory – Fabricate the Working Cast

- Verify that the OverdenSURE Abutment Analogs [ODS-AA01] are secure in the Impression Copings and pour the model.



Laboratory – Record Base and Occlusion Rim

- Fabricate a record base and wax occlusion rim and send to the clinician.



Clinician – Occlusion Record

- Place the occlusion rim in the mouth and take the interocclusal records.



Laboratory – Articulate

- Articulate casts using the interocclusal record.
- Set the denture teeth on the record base for the wax try-in.



Clinician – Denture Try-In

- Place the wax denture in the mouth. Verify occlusion, aesthetics, and phonetics. Make any necessary adjustments.
- If major adjustments are necessary, make a new interocclusal record and return it to the laboratory for a new articulator mounting, wax denture set-up, and try-in.



Laboratory – Process the Removable Prosthesis

- After the wax try-in is verified, flask the denture on the cast.
- Boil out the wax and separate the flask.
- Place the OverdenSURE White Spacers over Abutment Analogs to prevent acrylic resin from flowing under the Housing Caps.
- Place the OverdenSURE Housing Caps [ODS-RH01] with the pre-loaded Black Inserts on the Abutment Analogs and press down firmly to make sure they are fully seated.
- Process the denture following standard procedures.




Clinician – Remove OverdenSURE Black Inserts

- Using the OverdenSURE Removal and Insertion Tool [ODS-IRTOOL], remove the Black Inserts from each Housing Cap.
- To do so, make sure the Removal and Insertion Tool is set to "removal mode." This is achieved by turning the golden head of the tool two full rotations counter-clockwise. You will notice the blue tip recess into the tool.



- Insert the tool into the Black Insert and tilt it so the edge of the tip engages the inner lip of the insert.
- Once the tool has full engagement pull it out.
- To remove the Black Insert from the tool, turn the golden head of the tool two full rotations clockwise. The blue tip will emerge and push off the Black Insert.



Clinician – Place Final OverdenSURE Retention Inserts


- Ensure the Removal and Insertion Tool [ODS-IRTOOL] is set to "insertion mode." To do so, completely unscrew the golden head from the blue portion of the tool.
- Select the appropriate Retention Insert and place it on the insertion tip of the tool.
- Firmly insert each final Retention Insert into the empty Housing Caps, and verify they sit flush with the rim of the Housing Caps.

Note: In all geographies with the exception of the U.S., an OverdenSURE Standard Abutment Retention Kit [ODS-RIKIT2001] is included with every abutment.

One Blue Retention Insert (Extra Light), one Pink Retention

Insert (Light), one Clear Retention Insert (Medium), and one Red Retention Insert (Extra Light; Extended Range) are included within this kit and packaged with every abutment.

Note: In the U.S., OverdenSURE Abutments are packaged individually and retention inserts (ODSZ-RIKIT2001 OR ODSZ-RIKIT4001) must be ordered separately.



Deliver the Removable Prosthesis

- Seat the prosthesis on the OverdenSURE Abutments and press down so that the Retention Inserts snap onto the abutments.
- Verify the occlusion.
- Have the patient remove and insert the prosthesis to ensure the retention is adequate. If the retention level needs to be adjusted, select the appropriate Retention Insert and repeat the "Selecting the OverdenSURE Retention Inserts" steps.
- Instruct the patient on proper hygiene care and required maintenance visits.



Retention Kit ODSZ-RIKIT2001 not included, please order separately. Available only in the U.S.

Ordering Info

OverdenSURE Attachments for ZimVie Implant Systems (U.S. markets only)

TSX and Eztetic Implant Systems

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
NP (2.9 mmD)	ODSA-NP201	ODSA-NP202	ODSA-NP203	ODSA-NP204	ODSA-NP205	ODSA-NP206
3.5 mmD	ODSA-TSV301	ODSA-TSV302	ODSA-TSV303	ODSA-TSV304	ODSA-TSV305	ODSA-TSV306
4.5 mmD	ODSA-TSV401	ODSA-TSV402	ODSA-TSV403	ODSA-TSV404	ODSA-TSV405	ODSA-TSV406



TSV Implant Systems

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
3.5 mmD	ODSA-TSV301	ODSA-TSV302	ODSA-TSV303	ODSA-TSV304	ODSA-TSV305	ODSA-TSV306
4.5 mmD	ODSA-TSV401	ODSA-TSV402	ODSA-TSV403	ODSA-TSV404	ODSA-TSV405	ODSA-TSV406
5.7 mmD	ODSA-TSV501	ODSA-TSV502	ODSA-TSV503	ODSA-TSV504	ODSA-TSV505	ODSA-TSV506



Certain Implant Systems

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
3.4 mmD	ODSA-CERT301	ODSA-CERT302	ODSA-CERT303	ODSA-CERT304	ODSA-CERT305	ODSA-CERT306
4.1 mmD	ODSA-CERT401	ODSA-CERT402	ODSA-CERT403	ODSA-CERT404	ODSA-CERT405	ODSA-CERT406
5.0 mmD	ODSA-CERT501	ODSA-CERT502	ODSA-CERT503	ODSA-CERT504	ODSA-CERT505	ODSA-CERT506



External Hex Implant Systems

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
3.4 mmD	ODSA-EXHEX301	ODSA-EXHEX302	ODSA-EXHEX303	ODSA-EXHEX304	ODSA-EXHEX305	ODSA-EXHEX306
4.1 mmD	ODSA-EXHEX401	ODSA-EXHEX402	ODSA-EXHEX403	ODSA-EXHEX404	ODSA-EXHEX405	ODSA-EXHEX406
5.0 mmD	ODSA-EXHEX501	ODSA-EXHEX502	ODSA-EXHEX503	ODSA-EXHEX504	ODSA-EXHEX505	ODSA-EXHEX506



Retention Kit ODSZ-RIKIT2001 not included, please order separately. Available only in the U.S.

OverdenSURE Attachments for ZimVie Implant Systems (all markets outside the U.S.)

TSX and Eztetic Implant Systems

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
NP (2.9 mmD)	ODS-NP201	ODS-NP202	ODS-NP203	ODS-NP204	ODS-NP205	ODS-NP206
3.5 mmD	ODS-TSV301	ODS-TSV302	ODS-TSV303	ODS-TSV304	ODS-TSV305	ODS-TSV306
4.5 mmD	ODS-TSV401	ODS-TSV402	ODS-TSV403	ODS-TSV404	ODS-TSV405	ODS-TSV406



TSV Implant Systems

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
3.5 mmD	ODS-TSV301	ODS-TSV302	ODS-TSV303	ODS-TSV304	ODS-TSV305	ODS-TSV306
4.5 mmD	ODS-TSV401	ODS-TSV402	ODS-TSV403	ODS-TSV404	ODS-TSV405	ODS-TSV406
5.7 mmD	ODS-TSV501	ODS-TSV502	ODS-TSV503	ODS-TSV504	ODS-TSV505	ODS-TSV506



Certain Implant Systems

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
3.4 mmD	ODS-CERT301	ODS-CERT302	ODS-CERT303	ODS-CERT304	ODS-CERT305	ODS-CERT306
4.1 mmD	ODS-CERT401	ODS-CERT402	ODS-CERT403	ODS-CERT404	ODS-CERT405	ODS-CERT406
5.0 mmD	ODS-CERT501	ODS-CERT502	ODS-CERT503	ODS-CERT504	ODS-CERT505	ODS-CERT506



External Hex Implant Systems

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
3.4 mmD	ODS-EXHEX301	ODS-EXHEX302	ODS-EXHEX303	ODS-EXHEX304	ODS-EXHEX305	ODS-EXHEX306
4.1 mmD	ODS-EXHEX401	ODS-EXHEX402	ODS-EXHEX403	ODS-EXHEX404	ODS-EXHEX405	ODS-EXHEX406
5.0 mmD	ODS-EXHEX501	ODS-EXHEX502	ODS-EXHEX503	ODS-EXHEX504	ODS-EXHEX505	ODS-EXHEX506



Tapered SwissPlus Implant System

Seating Surface	Collar Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
3.8 mm	ODS-SWP301	ODS-SWP302	ODS-SWP303	ODS-SWP304	ODS-SWP305	ODS-SWP306
4.8 mm	ODS-SWP401	ODS-SWP402	ODS-SWP403	ODS-SWP404	ODS-SWP405	ODS-SWP406



Retention Kit ODS-RIKIT2001 included. Available in all active geographies except the U.S.

OverdenSURE Attachments for Multi-Platform Implant Systems (U.S. markets only)

DAEV-Conical compatible with Astra Tech Implant System™ EV			
Collar Height	Seating Surface		
	3.6 mm	4.2 mm	4.8 mm
1.0 mm	ODSA-DAEV361	ODSA-DAEV421	ODSA-DAEV481
2.0 mm	ODSA-DAEV362	ODSA-DAEV422	ODSA-DAEV482
3.0 mm	ODSA-DAEV363	ODSA-DAEV423	ODSA-DAEV483
4.0 mm	ODSA-DAEV364	ODSA-DAEV424	ODSA-DAEV484



NRT-Trilobe compatible with NobelReplace Select™		
Collar Height	Seating Surface	
	3.5 mm	4.3 mm
1.0 mm	ODSA-NRT351	ODSA-NRT431
2.0 mm	ODSA-NRT352	ODSA-NRT432
3.0 mm	ODSA-NRT353	ODSA-NRT433
4.0 mm	ODSA-NRT354	ODSA-NRT434



DATX-Internal compatible with Astra Tech Osseospeed™			
Collar Height	Seating Surface		
	3.0 mm	3.5/4.0 mm	4.5/5.0 mm
1.0 mm	-	ODSA-DATX351	ODSA-DATX451
2.0 mm	ODSA-DATX302	ODSA-DATX352	ODSA-DATX452
3.0 mm	ODSA-DATX303	ODSA-DATX353	ODSA-DATX453
4.0 mm	ODSA-DATX34	ODSA-DATX354	ODSA-DATX454



NAC-Conical compatible with NobelActive®/Replace®		
Collar Height	Seating Surface	
	3.5 mm	4.3/5.0 mm
1.0 mm	ODSA-NAC351	ODSA-NAC431
2.0 mm	ODSA-NAC352	ODSA-NAC432
3.0 mm	ODSA-NAC353	ODSA-NAC433
4.0 mm	ODSA-NAC354	ODSA-NAC434



BHI-Internal compatible with BioHorizons® Internal		
Collar Height	Seating Surface	
	3.5 mm	4.5 mm
1.0 mm	ODSA-BHI351	ODSA-BHI451
2.5 mm	ODSA-BHI3525	ODSA-BHI4525
3.5 mm	ODSA-BHI3535	ODSA-BHI4535
4.5 mm	ODSA-BHI3545	ODSA-BHI4545
5.5 mm	ODSA-BHI3555	-



NB-External compatible with Nobel Brånemark®		
Collar Height	Seating Surface	
	3.5 mm	4.1 mm
1.0 mm	ODSA-NB351	ODSA-NB411
2.0 mm	ODSA-NB352	ODSA-NB412
3.0 mm	ODSA-NB353	ODSA-NB413
4.0 mm	-	ODSA-NB414
5.0 mm	-	ODSA-NB415



CCA-Internal compatible with Camlog®		
Collar Height	Seating Surface	
	3.8 mm	4.3 mm
1.0 mm	ODSA-CCA381	ODSA-CCA431
2.0 mm	ODSA-CCA382	ODSA-CCA432
3.0 mm	ODSA-CCA383	ODSA-CCA433
4.0 mm	ODSA-CCA384	ODSA-CCA434



SBL-Conical compatible with Straumann® Bone level		
Collar Height	Seating Surface	
	3.3 mm	4.1/4.8 mm
1.0 mm	-	ODSA-SBL411
2.0 mm	ODSA-SBL332	ODSA-SBL412
3.0 mm	ODSA-SBL333	ODSA-SBL413
4.0 mm	ODSA-SBL334	ODSA-SBL414



DA-internal compatible with Dentsply Ankylos® C/X	
Collar Height	3.5/4.5 mm
2.0 mm	ODSA-DA352
3.0 mm	ODSA-DA353
4.0 mm	ODSA-DA354



STL-Tissue compatible with Straumann® Tissue level & synOcta	
Collar Height	4.1/4.8 mm
1.0 mm	ODSA-STL411
2.0 mm	ODSA-STL412
3.0 mm	ODSA-STL413
4.0 mm	ODSA-STL414



NG-Conical compatible with Neodent® GM	
Collar Height	Seating Surface
	3.0 mm
0.8 mm	ODSA-NG308
1.5 mm	ODSA-NG315
2.5 mm	ODSA-NG325
3.5 mm	ODSA-NG335
4.5 mm	ODSA-NG345
5.5 mm	ODSA-NG355




Retention Kit ODSZ-RIKIT2001 not included, please order separately. Available only in the U.S.

OverdenSURE Attachment Components

MSHI-Internal compatible with Mis® Seven		
Collar Height	Seating Surface	
	3.75/4.2 mm	5.0/6.0 mm
1.0 mm	ODSA-MSH341	ODSA-MSH451
2.5 mm	ODSA-MSH3425	ODSA-MSH4525
3.5 mm	ODSA-MSH3435	ODSA-MSH4535
4.5 mm	ODSA-MSH3445	ODSA-MSH4545
5.5 mm	ODSA-MSH3455	-



Description	Item No.
Retention Kit 20° Qty 1 Set	ODSZ-RIKIT2001
Retention Kit 20° Qty 2 Sets	ODSZ-RIKIT2002
Retention Kit 40° Qty 1 Set	ODSZ-RIKIT4001
Retention Kit 40° Qty 2 Sets	ODSZ-RIKIT4002
Blue Retention Insert 20° Qty 4	ODSZ-RIB2004
Pink Retention Insert 20° Qty 4	ODSZ-RIP2004
Clear Retention Insert 20° Qty 4	ODSZ-RIC2004
Red Retention Insert 40° Qty 4	ODSZ-RIR4004
Yellow Retention Insert 40° Qty 4	ODSZ-RIY4004
Green Retention Insert 40° Qty 4	ODSZ-RIG4004



OverdenSURE Accessories

Description	Item No.
Abutment Bar Attachment Qty 1	ODS-BARATT01
Abutment Bar Attachment Qty 2	ODS-BARATT02
Abutment Bar Attachment Qty 4	ODS-BARATT04
Analog Qty 1	ODS-AA01
Analog Qty 4	ODS-AA04
Analog Qty 10	ODS-AA010
Blockout Spacer Qty 20	ODS-BS020
Impression Coping Qty 1	ODS-IC01
Impression Coping Qty 4	ODS-IC04
Impression Coping Qty 10	ODS-IC010
Retention Housing Cap Qty 1	ODS-RH01
Retention Housing Cap Qty 4	ODS-RH04
Retention Housing Cap Qty 10	ODS-RH010
Scan Body Qty 10*	ODS-SB010



OverdenSURE Tools

Description	Item No.
Abutment Driver	ODS-DRVR
Abutment Driver for 7 mm (D) Straumann Torque Wrench	ODS-AD7ST
Abutment Driver for 8 mm (D) Noble Torque Wrench	ODS-AD8NB
ISO Latch, 4x4 Driver	ISO4X4A
Removal and Insertion Tool	ODS-IRTOOL



* Available in April 2025

OverdenSURE Attachments for Multi-Platform Implant Systems (all markets outside U.S.)

DAEV-Conical compatible with Astra Tech Implant System™ EV			
Collar Height	Seating Surface		
	3.6 mm	4.2 mm	4.8 mm
1.0 mm	ODS-DAEV361	ODS-DAEV421	ODS-DAEV481
2.0 mm	ODS-DAEV362	ODS-DAEV422	ODS-DAEV482
3.0 mm	ODS-DAEV363	ODS-DAEV423	ODS-DAEV483
4.0 mm	ODS-DAEV364	ODS-DAEV424	ODS-DAEV484



NRT-Trilobe compatible with NobelReplace Select™		
Collar Height	Seating Surface	
	3.5 mm	4.3 mm
1.0 mm	ODS-NRT351	ODS-NRT431
2.0 mm	ODS-NRT352	ODS-NRT432
3.0 mm	ODS-NRT353	ODS-NRT433
4.0 mm	ODS-NRT354	ODS-NRT434



DATX-Internal compatible with Astra Tech Osseospeed™			
Collar Height	Seating Surface		
	3.0 mm	3.5/4.0 mm	4.5/5.0 mm
1.0 mm	-	ODS-DATX351	ODS-DATX451
2.0 mm	ODS-DATX302	ODS-DATX352	ODS-DATX452
3.0 mm	ODS-DATX303	ODS-DATX353	ODS-DATX453
4.0 mm	ODS-DATX34	ODS-DATX354	ODS-DATX454



NAC-Conical compatible with NobelActive® /Replace®		
Collar Height	Seating Surface	
	3.5 mm	4.3/5.0 mm
1.0 mm	ODS-NAC351	ODS-NAC431
2.0 mm	ODS-NAC352	ODS-NAC432
3.0 mm	ODS-NAC353	ODS-NAC433
4.0 mm	ODS-NAC354	ODS-NAC434



BHI-Internal compatible with BioHorizons® Internal		
Collar Height	Seating Surface	
	3.5 mm	4.5 mm
1.0 mm	ODS-BHI351	ODS-BHI451
2.5 mm	ODS-BHI3525	ODS-BHI4525
3.5 mm	ODS-BHI3535	ODS-BHI4535
4.5 mm	ODS-BHI3545	ODS-BHI4545
5.5 mm	ODS-BHI3555	-



NB-External compatible with Nobel Brånemark®		
Collar Height	Seating Surface	
	3.5 mm	4.1 mm
1.0 mm	ODS-NB351	ODS-NB411
2.0 mm	ODS-NB352	ODS-NB412
3.0 mm	ODS-NB353	ODS-NB413
4.0 mm	ODS-NB354	ODS-NB414
5.0 mm	-	ODS-NB415



CCA-INTERNAL compatible with Camlog®		
Collar Height	Seating Surface	
	3.8 mm	4.3 mm
1.0 mm	ODS-CCA381	ODS-CCA431
2.0 mm	ODS-CCA382	ODS-CCA432
3.0 mm	ODS-CCA383	ODS-CCA433
4.0 mm	ODS-CCA384	ODS-CCA434



SBL-Conical compatible with Straumann Bone level		
Collar Height	Seating Surface	
	3.3 mm	4.1/4.8 mm
1.0 mm	-	ODS-SBL411
2.0 mm	ODS-SBL332	ODS-SBL412
3.0 mm	ODS-SBL333	ODS-SBL413
4.0 mm	ODS-SBL334	ODS-SBL414



DA-internal compatible with Dentsply Ankylos® C/X	
Collar Height	3.5/4.5 mm
2.0 mm	ODS-DA352
3.0 mm	ODS-DA353
4.0 mm	ODS-DA354



STL-Tissue compatible with Straumann® Tissue level & synOcta	
Collar Height	4.1/4.8 mm
1.0 mm	ODS-STL411
2.0 mm	ODS-STL412
3.0 mm	ODS-STL413
4.0 mm	ODS-STL414



NG-Conical compatible with Neodent® GM	
Collar Height	3.0 mm
0.8 mm	ODS-NG308
1.5 mm	ODS-NG315
2.5 mm	ODS-NG325
3.5 mm	ODS-NG335
4.5 mm	ODS-NG345
5.5 mm	ODS-NG355



Retention Kit ODS-RIKIT2001 included. Available in all active geographies except the U.S.

OverdenSURE Attachment Components

MSHI-Internal compatible with Mis® Seven		
Collar Height	Seating Surface	
	3.75/4.2 mm	5.0/6.0 mm
1.0 mm	ODS-MSH341	ODSA-MSH451
2.5 mm	ODS-MSH3425	ODSA-MSH4525
3.5 mm	ODS-MSH3435	ODSA-MSH4535
4.5 mm	ODS-MSH3445	ODSA-MSH4545
5.5 mm	ODS-MSH3455	-



Description	Item No.
Retention Kit 20° Qty 1 Set	ODS-RIKIT2001
Retention Kit 20° Qty 2 Sets	ODS-RIKIT2002
Retention Kit 40° Qty 1 Set	ODS-RIKIT4001
Retention Kit 40° Qty 2 Sets	ODS-RIKIT4002
Blue Retention Insert 20° Qty 4	ODS-RIB2004
Pink Retention Insert 20° Qty 4	ODS-RIP2004
Clear Retention Insert 20° Qty 4	ODS-RIC2004
Red Retention Insert 40° Qty 4	ODS-RIR4004
Yellow Retention Insert 40° Qty 4	ODS-RIY4004
Green Retention Insert 40° Qty 4	ODS-RIG4004



OverdenSURE Accessories

Description	Item No.
Abutment Bar Attachment Qty 1	ODS-BARATT01
Abutment Bar Attachment Qty 2	ODS-BARATT02
Abutment Bar Attachment Qty 4	ODS-BARATT04
Analog Qty 1	ODS-AA01
Analog Qty 4	ODS-AA04
Analog Qty 10	ODS-AA010
Blockout Spacer Qty 20	ODS-BS020
Impression Coping Qty 1	ODS-IC01
Impression Coping Qty 4	ODS-IC04
Impression Coping Qty 10	ODS-IC010
Retention Housing Cap Qty 1	ODS-RH01
Retention Housing Cap Qty 4	ODS-RH04
Retention Housing Cap Qty 10	ODS-RH010
Scan Body Qty 10*	ODS-SB010



OverdenSURE Tools

Description	Item No.
Abutment Driver	ODS-DRVR
Abutment Driver for 7 mm (D) Straumann Torque Wrench	ODS-AD7ST
Abutment Driver for 8 mm (D) Noble Torque Wrench	ODS-AD8NB
ISO Latch, 4x4 Driver	ISO4X4A
Removal and Insertion Tool	ODS-IRTOOL



* Available in April 2025

References:

1. Oh DJ, Kim HJ, Chung CH. A study on mechanical properties of TiN, ZrN and WC coated film on the titanium alloy surface. The Journal of Korean Academy of Prosthodontics. 2006;44(6):740-50.
2. Ul-Hamid A. Synthesis, microstructural characterization and nanoindentation of Zr, Zr-nitride and Zr-carbonitride coatings deposited using magnetron sputtering. J Adv Res. 2020 Nov 25;29:107-119.
3. Grössner-Schreiber B, Griepentrog M, Haustein I, Müller WD, Lange KP, Briedigkeit H, Göbel UB. Plaque formation on surface modified dental implants. An in vitro study. Clin Oral Implants Res. 2001 Dec;12(6):543-51.
4. Groessner-Schreiber B, Hannig M, Dück A, Griepentrog M, Wenderoth DF. Do different implant surfaces exposed in the oral cavity of humans show different biofilm compositions and activities?. European journal of oral sciences. 2004 Dec;112(6):516-22.
5. Grössner-Schreiber B, Teichmann J, Hannig M, Dörfer C, Wenderoth DF, Ott SJ. Modified implant surfaces show different biofilm compositions under in vivo conditions. Clin Oral Implants Res. 2009 Aug;20(8):817-26.
6. Brunello G, Brun P, Gardin C, Ferroni L, Bressan E, Meneghello R, Zavan B, Sivoletta S. Biocompatibility and antibacterial properties of zirconium nitride coating on titanium abutments: An in vitro study. PLoS One. 2018 Jun 26;13(6):e0199591.
7. Zhurakivska K, Ciacci N, Troiano G, Caponio VC, Scrascia R, Pallecchi L, Lo Muzio L, Arena F. Nitride-Coated and Anodic-Oxidized Titanium Promote a Higher Fibroblast and Reduced Streptococcus gordonii Proliferation Compared to the Uncoated Titanium. Prosthesis. 2020 Dec;2(4):333-9.

For more information, visit [ZimVie.com](https://www.zimvie.com)

ZimVie US Corp LLC
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272
dentalCS@ZimVie.com



ASTRA TECH™, ANKYLOS® C/X, ASTRATECH implant system™ EV are registered trademarks of DENTSPLY SIRONA Implants. BioHorizons® is a registered trademark of BioHorizons, Inc. Brånemark System®, NobelReplace® and Replace Select™, NobelActive™, NobelReplace® NobelSpeedy™ and Multi-Unit® are registered trademarks of Nobel Biocare AB. Straumann® and synOcta® are registered trademarks of Straumann Holding AG. , SynOcta(R), Neodent(R) and Grand Morse™. Locator® is registered trademark of Zest Dental Solutions®.

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. and its affiliates; and all products are manufactured by one or more of the dental subsidiaries of ZimVie US Corp LLC (Biomet 3i, LLC, etc.) and distributed by ZimVie US Corp LLC and its authorized marketing partners. OverdenSURE is manufactured by Terrats Medical S.L. EU Medical Device classification: Overdenture Abutments (Class IIb), Dental Implants (Class IIb); Overdenture Attachments, Components and Tools (Class I). The Intended purpose of the OverdenSURE Attachment System is to be attached to dental implants and act as support to attach the prosthesis to allow the improvement or restoration. The intended purpose of Dental Implants is to replace missing individual teeth, multiple teeth, or to restore edentulous dental arches. The products here shown are not eligible for reimbursement in France. For other countries, please check the local laws and regulation. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZVINST0073-GLBL REV F 11/25 ©2025 ZimVie. All rights reserved.



The NB number applicable to each device can be found on the product label, if applicable.